

PORTABLE RAMPS



Distributed in SE Wisconsin through:



877-917-2273



SUITCASE® Signature Series®



The SUITCASE® Signature Series® portable wheelchair ramp, with its single-fold design, sets up quickly and can be carried like a suitcase when not in use! It features a no-pinch, non-protruding full-length live hinge and extruded tread for superior traction that distinguishes it from its competition. Finished with durable, ergonomically designed handles, this ramp was redesigned to be simply the best in the industry. Top Lip Extension (TLE) compatible.

Available in 2' – 8' lengths with an 800-pound weight capacity.



Features a no-pinch, non-protruding live hinge for safe, silent, and maintenance free operation




Ergonomically designed flexible, non-breakable handles offer convenient carrying and a comfortable grip



Bottom transition plates self-adjust to accommodate uneven terrain

TOP LIP EXTENSION (TLE)

The TLE will increase the reach of the SUITCASE® Signature Series® by extending the top lip from 3" to 9"! The TLE is unique because it can be added or removed as needed, eliminating the need for a large ramp inventory. Shipped as a pair.



With the use of a TLE, the rear bumper is easily cleared



Patented

SUITCASE® Signature Series®

SPECIFICATIONS	SUITCASE SS2	SUITCASE SS3	SUITCASE SS4	SUITCASE SS5	SUITCASE SS6	SUITCASE SS7	SUITCASE SS8	TLE
Made of aluminum	•	•	•	•	•	•	•	•
Non-skid surface	•	•	•	•	•	•	•	•
Usable size	2' x 30"	3' x 30"	4' x 30"	5' x 30"	6' x 30"	7' x 30"	8' x 30"	8½" x 13⅞" x 2⅜"
Folded size (L x W x H)	2' x 15" x 4"	3' x 15" x 4"	4' x 15" x 4"	5' x 15" x 4"	6' x 15" x 4"	7' x 15" x 4"	8' x 15" x 4"	n/a
Built-in carrying handle	•	•	•	•	•	•	•	n/a
Weight (lbs.) *	9	14	18	25	30	36	42	5 per set
Weight capacity (lbs.)	800	800	800	800	800	800	800	**500
For wheelchairs & scooters	•	•	•	•	•	•	•	•
Usage guidelines	Curbs, 1 step	Curbs, 1 step	Curbs, 1-2 steps	Curbs, 1-2 steps, some minivans	1-3 steps, some minivans & SUVs	1-3 steps, some minivans & SUVs	2-3 steps, some minivans & SUVs	n/a

* Weights are approximate and may vary. Consult your equipment's owner guide for proper degree of incline. Never exceed its recommendations.

** Ramp weight capacity is reduced when used with the TLE.



Bariatric weight capacity = 800 lbs. or greater



UPS Shippable




SUITCASE® Advantage Series®



The Advantage Series® is a SUITCASE ramp you can separate, making it the easiest to carry and move. By simply removing two hinge pins, it quickly separates into two lightweight sections, each with its own ergonomically designed handle.

This ultra-strong, single-fold ramp provides the safety and strength required while remaining light and portable. Features an applied skid-resistant tape with self-adjusting bottom transition plates for easy conversion from ramp to ground.

Manufactured of welded aircraft-grade aluminum.




Bottom transition plates adjust to uneven terrain



Ergonomically designed flexible, non-breakable handles offer convenient carrying and a comfortable grip

TOP LIP EXTENSION (TLE)

Extend the usefulness of your SUITCASE® Advantage Series® portable wheelchair ramp with the Top Lip Extension (TLE). This optional plate increases the top lip from 3" to 9", allowing the ramp to clear the rear bumper of SUVs and vans. Shipped as a pair.



Patented

SUITCASE® Advantage Series®

SPECIFICATIONS	SUITCASE AS2	SUITCASE AS3	SUITCASE AS4	SUITCASE AS5	SUITCASE AS6	TLE
Made of aluminum	•	•	•	•	•	•
Non-skid surface	•	•	•	•	•	•
Usable size	2' x 29½"	3' x 29½"	4' x 29½"	5' x 29½"	6' x 29½"	8½" x 13⅞" x 2⅝"
Folded size (L x W x H)	2' x 15" x 4½"	3' x 15" x 4½"	4' x 15" x 4½"	5' x 15" x 4½"	6' x 15" x 4½"	n/a
Built-in carrying handle	•	•	•	•	•	n/a
Weight (lbs.) *	6 ea. (2 sections)	8 ea. (2 sections)	10 ea. (2 sections)	13 ea. (2 sections)	15 ea. (2 sections)	5 per set
Weight capacity (lbs.)	800	800	800	800	800	**500
For wheelchairs & scooters	•	•	•	•	•	•
Usage guidelines	Curbs, 1 step	Curbs, 1 step	Curbs, 1-2 steps	Curbs, 1-2 steps, some minivans	1-3 steps, some minivans & SUVs	n/a

* Weights are approximate and may vary.

Weights listed are per ramp section (each section can be carried, assembled, and stored separately).

Consult your equipment's owner guide for proper degree of incline. Never exceed its recommendations.

** Ramp weight capacity reduced when used with the TLE.



Bariatric weight capacity = 800 lbs. or greater



UPS Shippable

TRIFOLD® Advantage Series®



The TRIFOLD® Advantage Series®, with its unique 3-fold design, offers the length required for wheelchairs and scooters to easily access steps, vehicles, and raised landings. The ramp is designed to be used and carried as one unit or can be separated into two individual sections.

The TRIFOLD features interlocking brackets to provide smoother operation and reduced pinch points. It features an applied skid-resistant surface and self-adjusting bottom transition plates for easy conversion from ramp to ground.


Manufactured of welded aircraft-grade aluminum.



Bottom transition plates adjust to uneven terrain



Ergonomically designed handles offer convenient carrying and a comfortable grip



Separates for transport and storage



Hinge

TOP LIP EXTENSION (TLE)

Extend the usefulness of your TRIFOLD® Advantage Series® portable wheelchair ramp with the Top Lip Extension (TLE). This optional plate increases 3" top lip to reach 9", allowing the ramp to clear the rear bumper of SUVs and vans by simply attaching to the top lip. Shipped as a pair.



Patented

TRIFOLD® Advantage Series®

SPECIFICATIONS	TRIFOLD AS5	TRIFOLD AS6	TRIFOLD AS7	TRIFOLD AS8	TRIFOLD AS10	TLE
Made of aluminum	•	•	•	•	•	•
Non-skid surface	•	•	•	•	•	•
Usable size	5' x 29½"	6' x 29½"	7' x 29½"	8' x 29½"	10' x 29½"	8½" x 13 ⁷ / ₈ " x 2 ³ / ₈ "
Folded size (L x W x H)	30¾" x 15" x 11½"	36¾" x 15" x 11½"	42¾" x 15" x 11½"	48¾" x 15" x 11½"	60¾" x 15" x 11½"	n/a
Built-in carrying handle	•	•	•	•	•	n/a
Weight (lbs.) *	16 ea. (2 sections)	19 ea. (2 sections)	22 ea. (2 sections)	25 ea. (2 sections)	30 ea. (2 sections)	5 per set
Weight capacity (lbs.)	800	800	800	800	800	**500
For wheelchairs & scooters	•	•	•	•	•	•
Usage guidelines	Curbs, 1-2 steps, some minivans	1-3 steps, minivans, some vans & SUVs	1-3 steps, minivans, some vans & SUVs	2-3 steps, minivans, some vans & SUVs	2-4 steps, minivans, some vans & SUVs	n/a

* Weights are approximate and may vary. Weights listed are per ramp section (each section can be carried, assembled, and stored separately). Consult your equipment's owner guide for proper degree of incline. Never exceed its recommendations.

** Ramp weight capacity reduced when used with the TLE.



Bariatric weight capacity = 800 lbs. or greater



UPS Shippable



Modular Entry Ramp



Ideal for sliding doors (two ramps shown)

The Transitions™ Modular Entry Ramp is a lightweight yet durable free-standing ramp designed for doorways, sliding glass doors, and raised landings. The modular design, with its interchangeable components, easily accommodates thresholds from ¾" to 6" high.


Easily add and remove sections to accommodate multiple thresholds! Ideal for use with scooters, wheelchairs, and walkers. Made of aluminum and is can be used both inside and outside.



Modular design features an adjustable top flap

Modular Entry Ramp							
SPECIFICATIONS	THRESH 1	THRESH 1.5	THRESH 2	THRESH 3	THRESH 4	THRESH 5	THRESH 6
Made of aluminum	•	•	•	•	•	•	•
Non-skid surface	•	•	•	•	•	•	•
Usable size (H x L x W)	1" x 6½" x 34"	1½" x 9½" x 34"	2" x 12" x 34"	3" x 17" x 34"	4" x 22¼" x 34"	5" x 27½" x 34"	6" x 32¾" x 34"
Weight (lbs.) *	3	5	5½	8	11	15	20
Weight capacity (lbs.)	600	600	600	600	600	600	600
For wheelchairs & scooters	•	•	•	•	•	•	•
Usage guidelines	¾" - 1¼" rise	1¼" - 1¾" rise	1¾" - 2¼" rise	2¾" - 3¼" rise	3¾" - 4¼" rise	4¾" - 5¼" rise	5¾" - 6¼" rise

* Weights are approximate and may vary. Consult your equipment's owner guide for proper degree of incline. Never exceed its recommendations.



Angled Entry Ramp



The Transitions™ Angled Entry Ramp is a self-supporting solution with adjustable legs to give you approximately 2" of flexibility to adapt to your desired heights.

This adjustment feature is especially helpful for locations where the surface may be slightly uneven.

Ideal for single step rises and accommodates doorways swinging both inward and out. Standard width of 36". All aluminum design, and features the extruded non-skid tread for superior traction.

STANDARD/SLIDING DOORS: The Transitions™ line of Entry Ramps are designed for both standard and sliding door applications. Often times, the ground is slightly sloped, and the ability to have adjustable feet is necessary... There's no easier solution than the Angled Entry Ramp. Simply measure the rise, and choose the height adjustment you need in the following Specifications chart.

* Weights are approximate and may vary.
Consult your equipment's owner guide for proper degree of incline.
Never exceed its recommendations.



Angled Entry Ramp

SPECIFICATIONS	THSS 12	THSS 24
Made of aluminum	•	•
Non-skid surface	•	•
Usable size (H x L x W)	12" x 36" x 1 1/4"	24" x 36" x 1 1/4"
Weight (lbs.) *	7	14
Weight capacity (lbs.)	600	600
For wheelchairs & scooters	•	•
Usage guidelines	1 1/4" - 2 1/4" rise	1 7/8" - 4 1/4" rise

Angled Entry Mat

features Beveled Edges

Simply place the Transitions™ Angled Entry Mat in front of the door to provide a smooth ground-to-sill transition. This slip-resistant rubber ramp is designed to be used on any surface and can easily be trimmed or notched for a custom fit.

Made of 100% recycled materials. Ideal for indoor and outdoor use and can easily be moved to another location if needed. Entry Mat does not interfere with door tracks, making it ideal for use with sliding glass doors.

Usable width is 36" and available in 1 1/2" and 2 1/2" heights. (see chart on next page)






Modular Entry Mat

Simply place the Transitions™ Modular Entry Mat against the door threshold (inside or outside) to provide a smooth ground-to-sill transition. Standard width is 48" (for wider entrances, place two ramps side-by-side).

Can be easily trimmed for lower heights and widths, or combined to reach greater heights. Rubber mats are slip-resistant and designed for years of continued use. Made of recycled rubber.



*This THR 250 is shown on top of a THR RISER 225



SPECIFICATIONS	Modular Entry Mat		Angled Entry Mat	
	THR 250-1	THR RISER 225-1PR	THRBE 150-1	THRBE 250-1
Made of recycled rubber	•	•	•	•
Non-skid surface	•	•	•	•
Usable size (H x L x W)	2 1/2" x 24" x 48"	2 1/4" x 24" x 24" (2 each)	1 1/2" x 14" x 36"	2 1/2" x 25" x 36"
Overall size (H x L x W)	2 1/2" x 24" x 48"	2 1/4" x 24" x 24" (2 each)	1 1/2" x 14" x 40"	2 1/2" x 25" x 43"
Weight (lbs.) *	40	23 1/2 per riser (2 each)	12	40
Weight capacity (lbs.)	850	850	850	850
For wheelchairs & scooters	•	•	•	•
Usage guidelines	2 1/2" rise**	Call customer service for details**	1 1/2" rise	2 1/2" rise

* Weights are approximate and may vary.

** Optional components available to reach heights greater than 2 1/2". Call for details. Consult your equipment's owner guide for proper degree of incline. Never exceed its recommendations.



Bariatric weight capacity = 800 lbs. or greater




UPS Shippable



Bariatric width = 48" or greater

This versatile Transitions™ Modular Entry Mat also adapts to your vehicle. Many vehicles have a 1"- 3" drop from entrance to cargo floor. The Modular Entry Mat easily trims to the contours of vehicle filling the void, providing a smooth transition for mobility equipment access.





EZ-ACCESS[®]

Moving you forward[®]

Distributed in SE Wisconsin through:



877-917-2273